

## Krebs-Henseleit Bicarbonate Buffer

*Cell Culture Medium*

产品编号: CM015

**产品简介:** 早在 19 世纪, 研究人员就探索如何在体外培养动物组织, 平衡盐 (balanced salt solution) 溶液是组织培养较早使用的溶液。1985 年 Sydney Ringer 在此基础上制备了完全由无机盐组分的盐溶液成功维持了哺乳动物心肌组织在体外的搏动 (contractility)。随后 Tyrode 制备了低特异性的、可广泛用于哺乳动物细胞工作的盐溶液 (Tyrode salt solution), 该培养基在当时被广为接受作为母液用以稀释其它比如蛋白质组分。从那时起, 多种平衡盐溶液开发用以组织和细胞培养, 根据其作用, 目前常用的平衡盐溶液分为以下几类:

- 1) 为细胞生长和增殖提供水分和无机盐类以满足正常细胞代谢;
- 2) 为细胞生长提供适宜的酸碱环境, 一般 pH 范围在 7.2-7.6 之间;
- 3) 提供氨基酸等营养素满足细胞生长和增殖代谢之需;
- 4) 添加碳水化合物 (glucose) 及代谢中间产物满足细胞的能量需求;
- 5) 为细胞提供适宜的渗透压环境。

**产品类型:** 无菌过滤即用型液体培养基

**渗透压:** 275±15

**酸碱度:** 7.3+0.2

**储存条件:** 4-8°C

**包装规格:** 500mL

**参考文献:** Krebs HA and Henseleit K. Studies on urea formation in the animal organism. Hoppe-Seylers Z. Physiol. Chem. 210:33-66, 1932.

**产品组分:**

Catalogue number	CM013	CM014	CM015
Ingredients	g/L	g/L	g/L
D-Glucose	2.0	1.8	
Magnesium Sulfate [Anhydrous]	0.141	0.141	0.141
Potassium Phosphate Monobasic	0.16	0.16	0.16
Potassium Chloride	0.35	0.35	0.35
Sodium Chloride	0.1	0.1	0.1
Calcium chloride dihydrate	0.373	0.373	0.373
Sodium bicarbonate	2.1	2.1	2.1
HEPES		5.958	5.958
<b>Specification</b>			
pH (after buffer)	7.3+0.2	7.3+0.2	7.3+0.2
Osmolality (mOsm)	270±15	290±15	275±15

FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC AND THERAPEUTIC PROCEDURES



M&C Gene Technology • Phone: (010)8205-7786 • (010)8693-7385 • Fax: (010)8205-9875

E-mail: order@macgene.com • Tech Support: support@macgene.com • URL: http://www.macgene.com